



# Tinnitus Therapy

**Sleep Sound Therapy**  
**Tinnitus Relief Devices**  
**Pillow Speakers**



# Tinnitus Therapy

## Sound Oasis Pillow Speakers



Sleep while listening to your favourite sounds! Turn any pillow into a Sleep Therapy Pillow with these ultra-thin stereo speakers. Simply place under your existing pillow and enjoy your favourite music or sounds for optimal relaxation and comfort. With the option of built-in amplification and two high fidelity, ultra-thin stereo speakers, the amplified pillow speaker is perfect for surrounding one's head with music and sound for enhanced therapy or privacy. Ideal for help with sleep or managing tinnitus.



- PS2** Stereo Pillow Speaker
- PS3** Amplified Stereo Pillow Speakers

## Sound Pillow



The Sound Pillow® looks and feels just like a normal pillow, but features built-in stereo speakers so that the sound from a bedside sound generator can be discreet and intimate. The Sound Pillow® is fitted with a cable and 3.5mm jack plug for plugging into the Sound Oasis® sound generator, a personal stereo or clock radio. Available with or without volume control. (Not compatible with the Naturcare® Sounds Relaxation Ball). Sound Pillow®: 59cm x 39cm x 16cm. Travel Sound Pillow®: 39cm x 24cm x 10cm.



Please note: This product is non-returnable.

- ACC386B** Sound Pillow
- ACC386C** Travel Sound Pillow

## Sound Oasis Sound Cards

Sound Oasis® Sound Cards allow Sound Oasis® owners to expand their collection of available sounds. By simply inserting the sound card into slot located on the base of the Sound Oasis®, users can supplement their unit's built-in sound offerings with 12 additional digitally recorded sounds (compatible with Sound Oasis® V3 only). Supplied as either a single card or a pack of 5.

### Ear Therapy - 12 Tracks

Ocean surf, Rain, Stream, White noise 4k, White noise 6k, Full spectrum, Thera1 1k-10k, Thera2 1k-10k, Thera3 1k-10k, Thera 2k-8k, Thera 2.5k-8k, Thera 2.5k-5k, Thera .9k-3.2k

### Spa Retreat - 12 Tracks

Hawaiian surf, Gentle rain, Waterfall, Island stream, Asian garden, Summer night, Unwind, Reflection, Serenity, Easy days, Late night, Stillness

### Nature Journey - 12 Tracks

Thunderstorm, Crashing waves, Forest rain, Waterfall, Gentle surf, Ocean rain, Wind, Secluded lake, Caribbean surf, Jungle stream, Summer night, Amazon lagoon

### Sleep Sounds For Baby - 12 Tracks

White noise, Ocean surf, Gentle rain, Summer night, Dolphins, Prenatal, Heartbeat, Womb, Lullaby, Car ride, Lullaby 2, Woodlands



- ACC385-C8** Ear Therapy
- ACC385-C9** Spa Retreat
- ACC385-C10** Nature Journey
- ACC385-C11** Sleep Sounds For Baby

- ACC385-C8-5** 5 x Ear Therapy
- ACC385-C9-5** 5 x Spa Retreat
- ACC385-C10-5** 5 x Nature Journey
- ACC385-C11-5** 5 x Sleep Sounds For Baby

## Sound Oasis Mini Sound Machine

£%

The World's Smallest Sound Machine® will help you relax, fall asleep easier, manage tinnitus, and even improve concentration during the day. An unsurpassed collection of therapeutic sounds, selected to provide you with the best possible solution for your needs are supplied on a custom micro SD card 15 hour battery life.



**ACC385Q** Mini Sound Machine

## Sound Oasis Mini White Noise Machine

£%

Enjoy soothing white noise wherever you go! The latest innovation from Sound Oasis® gives you white noise in its most portable form yet. Up to 20 hours of run time per charge.



**ACC385P** Mini White Noise Machine

## Sound Oasis Tinnitus

£%

The tinnitus sound machine from Sound Oasis gives you the best selection of sounds developed for tinnitus relief plus the ability to customise the sound card to your preference.

Featuring world leading sound quality, sleep & nature sounds now play at higher quality than a CD to immerse you in relaxing audio. Supplied with a UK power adaptor, can also be powered by 4 AA batteries, available separately.



**ACC385N** Sound Oasis Tinnitus

## Sound Oasis Bluetooth

£%

Sleep better with any one of 10 built-in sounds or stream via Bluetooth™ from your mobile device or computer. Perfect for improving sleep, blocking noise, and managing tinnitus. Features high quality speaker for clear sound quality with powerful bass. Powered by built-in rechargeable battery for wireless freedom or from USB. 30'/10m Bluetooth™ wireless range.



**ACC385L** Sound Oasis Bluetooth

## Sound Oasis Travel

£%

The Sound Oasis® Travel bedside desktop sound generator & alarm clock features 19 authentic nature sounds to soothe and relieve the symptoms of tinnitus at night.

Features a backlit, 12/24 hour alarm clock with a digital calendar, 4 position timer & memory button, as well as a headphone jack. Powered by AC adaptor (included). Travel pouch & time zone chart included. Available as a bundle with pillow speaker (PS2)



**ACC385H** Sound Oasis Travel  
**ACC385H-BUN** Sound Oasis Travel Bundle

# Tinnitus Therapy

## GTS-2000 DELUXE Glo to Sleep Therapy Mask



DELUXE Glo to Sleep is a safe, natural and effective way to help you get the sleep you need. Ancient wisdom joins with modern technology creating a revolutionary new way for you to clear your mind, quickly relax and get to sleep naturally. Sound Oasis® DELUXE Glo to Sleep will help you focus your mind on the process of relaxation and sleep - naturally and effectively. Features adjustable brightness with 3 auto dim options: ramp down in 10, 20 or 30 minutes.



**ACC385J** GTS-2000 DELUXE

## Naturcare Sounds Relaxation Ball



Aids tinnitus sufferers in getting a restful night's sleep by playing digital recordings of real, natural sounds. Choice of 7 digitally recorded sounds: Heartbeat, Sea Waves, Brook, Birds, Woods, White Noise, Rain.

Features a volume control, a 30 minute timer, as well as no tapes or moving parts to wear out. 4 x AA batteries included, with UK and EU power adaptors available separately. Bundle pack with UK adaptor also available.



**ACC323** Naturcare Sounds Relaxation Ball  
**ACC323A** Naturecare Bundle Pack

**ACC323B** UK Mains Adaptor  
**ACC323C** EU Mains Adaptor

## MM Series Tinnitus Relief Devices



MM1



MM2



MM6



MM8+



MM10

### MM1 Generic Fit Sound Generator

Straight stepped ITC white noise generator, with vent, volume control & type 10 battery (50dB also features output trimmer).

- MM1-WHITE** MM1 - 50dB Output
- MM1-BLUE** MM1 - 60dB Output
- MM1-RED** MM1 - 80dB Output
- MM1-GREEN** MM1 - 95dB Output

### MM2 Generic Fit Sound Generator

Straight contoured ITC white noise generator, with vent, volume control & type 10 battery (50dB also features output trimmer).

- MM2-WHITE** MM2 - 50dB Output
- MM2-BLUE** MM2 - 60dB Output
- MM2-RED** MM2 - 80dB Output
- MM2-GREEN** MM2 - 95dB Output

### MM6 Generic Fit Sound Generator

Helix fit white noise generator, with volume control & type 10 battery (50dB also features output trimmer). Available in left or right models. Also available as a DAC digitally programmable model.

Replacement tubing assembly also available supplied individually or in packs of 20.

- MM6-WHITE-L** 50dB Output - Left
- MM6-BLUE-L** 60dB Output - Left
- MM6-RED-L** 80dB Output - Left
- MM6-DACTRT-L** DAC TRT - Left
- MM6-WHITE-R** 50dB Output - Right
- MM6-BLUE-R** 60dB Output - Right
- MM6-RED-R** 80dB Output - Right
- MM6-DACTRT-R** DAC TRT - Right
- MM6TUBE-ASSY** Tubing - Single
- MM6TUBE-ASSY-20** Tubing - 20 Pack

### MM8+ Generic Fit Sound Generator

Generic fit for left or right ear, tapered shell, ITC (bullet style) white noise generator. Features volume control and vent. Available in 60dB output.

- MM8+** MM8+ - 60dB Output

### MM10 Generic Fit Sound Generator

Ultra-miniature straight contoured ITC white noise generator with volume control, vent & type 10 battery (50dB also features output trimmer). Also available as a DAC digitally programmable model.

- MM10-WHITE** MM10 - 50dB Output
- MM10-BLUE** MM10 - 60dB Output
- MM10-RED** MM10 - 80dB Output
- MM10-GREEN** 95dB Output
- MM10-DACTRT** DAC TRT

### Anti-Feedback Sleeve

Ideal for CIC hearing systems, for eliminating feedback in the ear. Made of thin, soft silicone which can be cut down to size as needed. Also used with themm range of tinnitus devices. Supplied individually or in packs of 50.



- ACC168** Single
- ACC168-50** 50 Pack

# Tinnitus Therapy

## Sueno Pro S Generic Tinnitus Therapy Device - Graphite

The Sueno Pro S is a compact & stylish tinnitus therapy device, powered by a 312 Zinc Air type battery with an average battery life of 120 hours. Featuring a maximum output of 70dB, a frequency range of 600 Hz–6400 Hz and three digital programs with adjustable noise in nine frequency channels, as well as three types of modulation (slow, medium and fast) offering a best-in-class solution for tinnitus therapies.

The Pro S features an easy to use rocker switch and nanoShield technology for enhanced water repellence. The device is supplied as a generic fit for either left or right ear and includes tubing and domes to suit all canal sizes.

**3 included programs:**  
Noise Only  
Pink Noise  
Waves



**SUENO-S-GRA** Sueno Pro S

## Sueno Pro R Generic Tinnitus Therapy Device - Graphite

The Sueno Pro R is a compact & stylish tinnitus therapy device, powered by a 312 Zinc Air type battery with an average battery life of 130 hours. Featuring a maximum output of 70dB, a frequency range of 100 Hz–8000 Hz and three digital programs with adjustable noise in nine frequency channels, as well as three types of modulation (slow, medium and fast) offering a best-in-class solution for tinnitus therapies.

The Pro R features an easy to use rocker switch and nanoShield technology for enhanced water repellence. The device is supplied as a left or right ear fit and includes receivers in 2 sizes, and domes to suit all canal sizes.

**3 included programs:**  
Noise Only  
Pink Noise  
Waves



Sueno Pro R Left

**SUENO-R-GRA-L** Sueno Pro R - Left  
**SUENO-R-GRA-R** Sueno Pro R - Right

## Sueno Pro S Mini Open Fit Kit

Mini Open Fit Kit For Sueno Pro S. Contains everything you need for the perfect fitting device.

Includes Open fitting domes 5mm, 7mm & 10mm (2 of each size), open tube cleaning wire (2 Pack), left open fitting tubes in sizes 1B & 2B, and right open fitting tubes in sizes 1B & 2B.



**PT-MOF-KIT** Mini Open Fit Kit

## Sueno Pro S Adaptors & Hooks

Open tube adaptors and hooks for Sueno Pro S.



**ACC505A** Open Tube Adaptor  
**ACC506A** Hook

## Sueno Pro R Receivers

Supplied individually.



**REC-S0L** Left - Size 0  
**REC-S1L** Left - Size 1  
**REC-S2L** Left - Size 2  
**REC-S3L** Left - Size 3

**REC-S0R** Right - Size 0  
**REC-S1R** Right - Size 1  
**REC-S2R** Right - Size 2  
**REC-S3R** Right - Size 3